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SAFETY REQUIREMENTS

The present manual is an integral part of the appliance and must be carefully read before using it as it gives important indications with regard to safety of installation, use and maintenance. Therefore preserve with care.

- Make sure the appliance is in good condition after you have removed the packing. In case of doubt do not use it and contact an assistance centre authorised by TERMOZETA.
- Do not leave the packing elements (plastic bags, foam polystyrene, nails, etc.) within reach of children or persons with limited mental faculties as they are potential sources of danger.
- This appliance was exclusively designed for domestic use to weigh people. Any other use is considered improper.
- TERMOZETA cannot be held responsible for any damage caused by improper, incorrect and negligent use.
- Contact a service centre authorised by the company for any repairs and ask for original spare parts.
- If you decide to no longer use the appliance, make any dangerous parts harmless, especially for children who may use the appliance to play with.
- Abide by the existing standards in your country for instructions on its disposal.
- Dear Customer, congratulations on your purchase of the Body Control 5000, bathroom scales of an extremely modern design.
- Body Control 5000 is characterised by an innovative technology that instantly measures your weight, the body fat and percentage of water contained in your body. The special measuring sensors combined with personal factors like sex, age and height make it possible for the Body Control 5000 bathroom scales to help you maintain your ideal weight.

INSTRUCTIONS FOR USE

ATTENTION

People with pace-makers or other electrical medical devices should not use these scales because they may interfere with their correct functioning.

NOTE

The following people should not use this appliance:

- people on dialysis
 - femmes enceintes
 - athletes with an intense training schedule
 - persons younger than 6 or older than 100
- The values measured by this appliance are indicative and do not have any medical value. Consult your doctor before embarking on any diet.

Installing the battery

The battery is not supplied with the scales. Use a type CR2032 3V battery. Position the scales with the rear facing upwards so as to be able to reach the battery housing. Remove the closure plate (Fig. 3). Insert the battery observing the polarity.

NOTE

If you do not plan on using the appliance for a long period of time, remove the battery.

Positioning

Position the appliance in a fresh dry area, avoid humid and hot places. Position the scales on a hard, flat surface without vibrations in order to ensure a safe and accurate measurement. Position the scales in an area not exposed to direct sunlight, having excessive humidity or with very sudden changes of temperature. Position the scales away from appliances that generate heat.

Features

The Body Control 5000 scales measure weight (kg, lb, st:lb), the percentage of fat mass (% fat), body hydration (% TBW), muscles (% **M**), bone weight (kg **ª**), the calorific value (KCAL) and the body mass index (BMI). The personal data (sex, age and height) can be pre-set and stored in specific profiles (P0, P1, etc.). The automatic shut-down function, during the programming, is enabled within 8 seconds if you do not touch any of the keys or buttons. After completing the measuring process, the scales show the results twice and then they shut off automatically. The scales also have a function that reports when the battery goes dead and when the allowed maximum weight (180 kg) is exceeded.

Settings of the scales

Before using the scales you have to set the desired unit of measure (pounds - lb, stones - st, kilograms - kg). Press the SET button. The scales turn on. Wait for the value "0.0" to appear. Press UP/DOWN to set the desired unit of measure.

To also view the percentage of fat mass and the other values, your data (sex, age and height) must be entered into one of the 10 profiles available, from P0 to P9, following the instructions below:

- Press the SET button. The first parameter to set will start flashing: P0 profile (Fig. 4).
- Press the UP/DOWN buttons to select the desired memory position and confirm by pressing the SET button. The second parameter to select will immediately start flashing on the display: sex (male **M** / female **F**) (Fig. 5).
- Set the sex (male **M** / female **F**) using the UP/DOWN buttons and confirm selection by pressing the SET button. Similar to what was described above, the next parameter to select will flash on the display: height (cm / ft/in) (Fig. 6).
- Set the height using the UP/DOWN buttons and confirm selection by pressing the SET button. The last parameter to set will immediately start flashing on the display: age. The age must fall between 6 and 100 years (Fig. 7).
- Set the age using the UP/DOWN buttons and confirm selection by pressing the SET button.

The procedure for setting the first profile that can be memorised is now complete.

It is possible to store up to 10 profiles, or change the parameters whenever necessary, by repeating all of the steps described above.

NOTE

If you make a mistake during the programming, repeat the steps from the selection of the personal profile.

Weight and body fat measuring

- The scales can be used like normal scales.
- Position the scales on a flat stable surface.
- Step onto the scales while still turned off and wait for your weight to be measured. The measured value flashes on the display. The scales automatically turn off after you step off of them.
- To also view the percentage of fat mass and the other values after having set all of your personal data, proceed as follows:
- Press the SET button and call up the set profile (from P0 to P9) using the UP/DOWN buttons.
- The display will show the data regarding the selected profile and "0.0" will appear immediately afterwards (Fig. 8). The scales are now ready to be used.

ATTENTION

Do not step onto the scales when they are damp or if your feet are wet or damp as you could slip. To avoid possible falls or injuries, do not step onto the scales plate edge.

NOTE

To ensure accurate measurement, step onto the scales with bare feet and make sure you correctly place your feet in contact with the metal strip of the sensors (2). Make sure your feet do not touch each other. Only in this way the sensors will function correctly and the scales will provide the data relating to body fat percentage, body water percentage and the following values. If you step onto the scales before "0.0" appears, the scales do not take the measurement. If you do not step onto the scales within about 8 seconds from when "0.0" appears, the scales automatically turn off.

The weight value appears in the centre of the display. Stay on the scales to get

the other values. The indicator (o) will scroll from left to right on the display to show that the scales are calculating the other values. At this point the percentage of fat mass (% fat), body hydration (% TBW), muscles (% **M**), bone weight (kg **ª**), the calorific value (KCAL) and the body mass index (BMI) will appear in succession on the display.

NOTE

Do not step down from the scales until the fat mass (% fat) and body hydration (% TBW) have been displayed. The indication of the weight and of the other detected values are repeated twice. Then the appliance turns off automatically.

Read the evaluation of the percentage data measured by the scales in the following tables while also referring to age and sex.

Age	Man			
	underweight	medium	overweight	obese
6-30	<16.4	16.5 - 22.0	22.1 - 30.4	>30.5
31-60	<17.0	17.1 - 22.4	22.5 - 31.4	>31.5
61-80	<17.6	17.7 - 23.0	23.1 - 32	>32.1

Age	Woman			
	underweight	medium	overweight	obese
6-30	<13.6	13.7 - 18.0	18.1 - 26.4	>26.5
31-60	3.0 - 14.0	14.1 - 18.6	18.7 - 27	>27.1
61-80	3.0 - 14.4	14.5 - 19.0	19.1 - 28	>28.1

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Remember that the measuring of the weight and the determination of the other values are based on standardised and calculated measuring methods with formulas perfected by in-depth scientific studies.

Habits and lifestyle affect the water levels, as well as the following factors:

- After a bath, the measures can be lower.
- After a meal or a drink, the measures can be higher.
- During the menstrual periods, the measures can fluctuate or be less precise.
- In case of illness or during a period of great physical activity, there can be a loss of water due to perspiration and dehydration.

Battery replacement

When the battery is dead, "LO" appears on the display. Immediately replace the battery to ensure an accurate measurement. If the scales do not turn on during use, it is necessary to replace the battery by removing it from its housing and putting in a new one, as described in the paragraph on battery installation.

NOTE

Removing the battery does not delete the stored settings.

Error messages

The scales may display some error messages:

Message on display	Meaning
LO	Low batteries
ERR	Weight exceeds 180 kg
Err2	The body fat and body water percentages cannot be calculated (if the body fat percentage is less than 5% or higher than 75%, the reading cannot be accurate). The measuring sensors are not working correctly. Try to position yourself correctly and if the error persists, clean the sensors.

Remove the battery for a few moments and put it back into place if there is any other error message or malfunctioning of the display.

CLEANING AND MAINTENANCE

Clean the scales with a soft dry cloth.

Do not attempt to disassemble the scales plate. Use alcohol in order to clean the sensors and a glass cleaner (applied to a cloth first) in order to keep them shiny. Soaps are not allowed.

ATTENTION

- Never use acid solutions or abrasives.
- Never immerse your scales in water or other liquids.

INFORMATION FOR USERS

In accordance with European Directives 2002/95/CE, 2002/96/EC and 2003/108/CE on the restriction of the use of dangerous substances in electric and electronic equipment as well as their waste disposal

The barred symbol of the rubbish bin shown on the equipment indicates that, at the end of its useful life, the product must be collected separately from other waste.

The adequate separate collection for the subsequent start-up of the equipment sent to be recycled, treated and disposal of in an environmentally compatible way contributes to preventing possible negative effects on the environment and health and optimises the recycling and reuse of components making up the apparatus.

Therefore, any products that have reached the end of their useful life must be given to waste disposal centres specialising in separate collection of waste electrical and electronic equipment, or given back to the retailer at the time of purchasing new similar equipment, on a one for one basis.

Abusive disposal of the product by the user involves application of the administrative sanctions according to the laws in force.

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